

F. Fuest¹, M. Schütte¹, T. Li², J. Pareja², A. Dreizler², B. Böhm²

¹LaVision GmbH, Anna-Vandenhoeck-Ring 19, D-37081 Goettingen, Germany ²TU Darmstadt, Germany

ffuest@lavisoin.de

Motivation

- Turbulent flow phenomena are three-dimensional in nature
- Planar OH-LIF is a common tool to investigate turbulent flame characteristics
- Information in 3rd dimension, however, is lost in planar techniques
- Turbulent flame features such as flame holes can only be characterized by fully three-dimensional measurements
- Tomographic OH-LIF imaging as an approach to yield the full three-dimensional OH distribution

Experimental setup




Figure 1: Experimental setup for tomographic LIF imaging.


Detection

- Four intensified CCD cameras in one plane at 45° angle separation
- UV lenses f/8 combined with LaVision Image doubler
- 8 simultaneous views

Excitation

- Frequency-doubled output of a dye laser tuned to excite the Q1(8) transition ($\lambda = 283.55 \text{ nm}$) of the $A^2\Sigma^+ \leftarrow X^2\Pi$ ($v'=1, v''=0$) band of hydroxyl radicals
- 20 mJ at probe volume of $3 \times 3 \times 3 \text{ cm}^3$

Flame and setup photographs



Spatial resolution and measurement uncertainty

Structural analysis and distribution of noise

- Residual noise in xy-planes of reconstructed single-shot (Fig. 2a, 2x2 binning) is obtained from subtraction of the same data but using a 9x9 Gaussian filter before reconstruction
- Structure of the noise pattern is analyzed by spatial autocorrelation in x- and y-direction (Fig. 2b) within the red rectangle in Fig. 2a
→ Spatial resolution limit due to noise is found at near-zero spatial correlation at 4-7 voxel = $450^3\text{-}1050^3 \mu\text{m}^3$
- Signal-to-noise ratio of SNR = 30 is determined from average count number of 5 and a standard deviation of 0.17 for the near Gaussian distributed noise within the red rectangle (Fig. 2c)




Figure 2: a) Residual spatial noise pattern in tomographic reconstruction of laminar Bunsen flame. b) Averaged autocorrelation function in x and y direction within the red rectangle. c) Distribution of residual spatial noise.

Impact of number of views on reconstruction in laminar Bunsen flame

- Integrated volume projections of shot-avg of reconstructed volume:



→ Reconstruction with only 2 views possible but notably distorted

Tomographic reconstruction

- Simultaneous Multiplicative Algebraic Reconstruction Technique (SMART)
- 100 iterations
- Computational time for 8 views @ 16 cores (3.10 GHz, 128 GB RAM):
 - 5 min for 100M voxel of $75^3 \mu\text{m}^3$, no binning,
 - 45 s for 12M voxel of $15^3 \mu\text{m}^3$, 2x2 binning





Figure 3: Reconstructed z-planes of turbulent Bunsen flame

Figure 4: 3D-Isosurfaces

Comparison with 2D PLIF measurements

- Tomographic LIF and PLIF at same location, see Fig. 5
- Laminar premixed methane flame, $d = 13 \text{ mm}$
- Center plane
- One side of OH-LIF profile
- 4 mm above the nozzle exit

Figure 5: Reconstructed single-shot laminar flame and location of extracted intensity profiles.



- Different filter sizes for PLIF and different volume-sheet thicknesses for tomographic LIF




Figure 6: Left: Shot-averaged profiles of TLIF compared to PLIF with different filter sizes. Right: Different TLIF sheet thickness. Distance between 50% intensity and jet centerline $\Delta x \approx 4.5 \text{ mm}$.

- PLIF using 5x5 filter shows similar gradient at large OH intensity
- Steeper gradient with decrease of sheet thickness

- Spatial width and fluctuation of intensity profiles




Figure 7: Left: Single-shot and average profile from PLIF. Right: TLIF. Length L of 50% intensity and its standard deviation σ .

- Slightly broader than PLIF but high resolution preserved

Turbulent lifted jet flame

- A lifted turbulent methane jet flame, jet diameter = 8mm, lift-off height $\approx 28 \text{ mm}$
- Central jet Reynolds number 5000, outlet velocity of 10 m/s
- Jet surrounded by a laminar air co-flow of 0.2 m/s
- 100 iterations, voxel of $75^3 \mu\text{m}^3$, no binning




Figure 8: Single-shot visualization of a jet flame by different views using Isosurfaces. Top-left: outer contour; top-right: section through center line; down: parallel horizontal and vertical cross sections.